

For Research Use Only

# PCSK6 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60595-1-PBS



## Basic Information

Catalog Number:

60595-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26288

GenBank Accession Number:

NM\_001291309

GeneID (NCBI):

5046

UNIPROT ID:

P29122

Full Name:

proprotein convertase  
subtilisin/kexin type 6

Calculated MW:

106 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G6G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50847-1, PCSK6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60595-1-PBS. Detection antibody: 60595-2-PBS. Standard:Ag26288. Range: 0.195-50 ng/mL



Cytometric bead array standard curve of MP50847-2, PCSK6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60595-3-PBS. Detection antibody: 60595-1-PBS. Standard:Ag26288. Range: 0.391-100 ng/mL