

For Research Use Only

TLR3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60598-1-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 60598-1-PBS | GenBank Accession Number: BC094737 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 7098 | CloneNo.: 2A1A7 |
| Source: Mouse | UNIPROT ID: O15455 | |
| Isotype: IgG1 | Full Name: toll-like receptor 3 | |
| Immunogen Catalog Number: AG12163 | Calculated MW: 904 aa, 104 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50854-1, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-2-PBS. Standard:Ag12163. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50854-2, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-3-PBS. Standard:Ag12163. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50854-3, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-4-PBS. Standard:Ag12163. Range: 0.781-50 ng/mL.