

For Research Use Only

TRIM13 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60599-1-PBS



Basic Information

Catalog Number: 60599-1-PBS	GenBank Accession Number: BC003579	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10206	CloneNo.: 1E9C1
Source: Mouse	UNIPROT ID: O60858	
Isotype: IgG2b	Full Name: tripartite motif-containing 13	
Immunogen Catalog Number: AG16805	Calculated MW: 407 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50856-1, TRIM13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60599-1-PBS. Detection antibody: 60599-2-PBS. Standard: Ag16805. Range: 3.125-50 ng/mL