

For Research Use Only

RPL28 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60607-2-PBS



Basic Information

Catalog Number: 60607-2-PBS	GenBank Accession Number: BC010173	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6158	CloneNo.: 1E11D9
Source: Mouse	UNIPROT ID: P46779	
Isotype: IgG1	Full Name: ribosomal protein L28	
Immunogen Catalog Number: AG9229	Calculated MW: 137 aa, 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50863-1, RPL28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60607-1-PBS. Detection antibody: 60607-2-PBS. Standard: Ag9229. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50863-2, RPL28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60607-2-PBS. Detection antibody: 60607-3-PBS. Standard: Ag9229. Range: 0.195-100 ng/mL