

For Research Use Only

# PPP2R2A Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60608-4-PBS



## Basic Information

Catalog Number:

60608-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9851

GenBank Accession Number:

BC041071

GeneID (NCBI):

5520

UNIPROT ID:

P63151

Full Name:

protein phosphatase 2 (formerly 2A),  
regulatory subunit B, alpha isoform

Calculated MW:

477 aa, 52 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C11C11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50864-3, PPP2R2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60608-3-PBS. Detection antibody: 60608-4-PBS. Standard: Ag9851. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50864-4, PPP2R2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60608-4-PBS. Detection antibody: 60608-5-PBS. Standard: Ag9851. Range: 0.781-100 ng/mL.