

For Research Use Only

CDO1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60610-4-PBS



Basic Information

Catalog Number: 60610-4-PBS	GenBank Accession Number: BC024241	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1036	CloneNo.: 2B9A6
Source: Mouse	UNIPROT ID: Q16878	
Isotype: IgG1	Full Name: cysteine dioxygenase, type I	
Immunogen Catalog Number: AG3294	Calculated MW: 200 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50867-3, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-1-PBS. Detection antibody: 60610-4-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL.