

For Research Use Only

RNF133 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60611-5-PBS



Basic Information

Catalog Number: 60611-5-PBS	GenBank Accession Number: BC022038	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 168433	CloneNo.: 2E4C10
Source: Mouse	UNIPROT ID: Q8WVZ7	
Isotype: IgG1	Full Name: ring finger protein 133	
Immunogen Catalog Number: AG14534	Calculated MW: 376 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50868-3, RNF133 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60611-5-PBS. Detection antibody: 60611-6-PBS. Standard: Ag14534. Range: 0.098-100 ng/mL.