

For Research Use Only

lama3 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number:60613-1-PBS



Basic Information

Catalog Number:

60613-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30144

GenBank Accession Number:

NM_198129

GeneID (NCBI):

3909

UNIPROT ID:

Q16787

Full Name:

laminin, alpha 3

Calculated MW:

367 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A8C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50873-1, lama3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60613-1-PBS. Detection antibody: 60613-2-PBS. Standard:Ag30144. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50873-2, lama3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60613-3-PBS. Detection antibody: 60613-1-PBS. Standard:Ag30144. Range: 0.391-100 ng/mL