

For Research Use Only

KIAA0776 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60614-1-PBS



Basic Information

Catalog Number:

60614-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23404

GenBank Accession Number:

BC036379

GeneID (NCBI):

23376

UNIPROT ID:

O94874

Full Name:

KIAA0776

Purification Method:

Protein G Magarose purification

CloneNo.:

3E9G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50874-1, KIAA0776 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60614-1-PBS. Detection antibody: 60614-2-PBS. Standard: Ag23404. Range: 0.195-200 ng/mL