

For Research Use Only

# PIP4K2C Monoclonal antibody, PBS Only (Detector)

Catalog Number:60615-2-PBS



## Basic Information

Catalog Number:

60615-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10989

GenBank Accession Number:

BC028596

GeneID (NCBI):

79837

UNIPROT ID:

Q8TBX8

Full Name:

phosphatidylinositol-5-phosphate 4-kinase, type II, gamma

Calculated MW:

421 aa, 47 kDa

Purification Method:

Protein G purification

CloneNo.:

2B6D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50878-1, PIP4K2C Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60615-1-PBS. Detection antibody: 60615-2-PBS. Standard: Ag10989. Range: 0.195-100 ng/mL.