

For Research Use Only

PIP4K2C Monoclonal antibody, PBS Only (Capture)

Catalog Number:60616-1-PBS



Basic Information

Catalog Number:

60616-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG11198

GenBank Accession Number:

BC028596

GeneID (NCBI):

79837

UNIPROT ID:

Q8TBX8

Full Name:

phosphatidylinositol-5-phosphate 4-kinase, type II, gamma

Calculated MW:

421 aa, 47 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3B11E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50878-2, PIP4K2C Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60616-1-PBS. Detection antibody: 60615-4-PBS. Standard: Ag10989. Range: 0.195-100 ng/mL.