

For Research Use Only

BFAR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60618-4-PBS



Basic Information

Catalog Number: 60618-4-PBS	GenBank Accession Number: BC003054	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51283	CloneNo.: 1D9D11
Source: Mouse	UNIPROT ID: Q9NZS9	
Isotype: IgG1	Full Name: bifunctional apoptosis regulator	
Immunogen Catalog Number: AG7755	Calculated MW: 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50881-3, BFAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60618-4-PBS. Detection antibody: 60618-2-PBS. Standard: Ag7755. Range: 0.098-25 ng/mL