

For Research Use Only

LRP5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60619-1-PBS



Basic Information

Catalog Number: 60619-1-PBS	GenBank Accession Number: BC150595	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4041	CloneNo.: 3B10B3
Source: Mouse	UNIPROT ID: O75197	
Isotype: IgG1	Full Name: low density lipoprotein receptor-related protein 5	
Immunogen Catalog Number: AG20043	Calculated MW: 1615 aa, 179 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50882-1, LRP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60619-1-PBS. Detection antibody: 60619-2-PBS. Standard: Ag20043. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50882-2, LRP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60619-1-PBS. Detection antibody: 60619-3-PBS. Standard: Ag20043. Range: 1.563-100 ng/mL.