

For Research Use Only

LRP5 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60619-5-PBS



Basic Information

Catalog Number: 60619-5-PBS	GenBank Accession Number: BC150595	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4041	CloneNo.: 2C12H5
Source: Mouse	UNIPROT ID: O75197	
Isotype: IgG2b	Full Name: Low density lipoprotein receptor-related protein 5	
Immunogen Catalog Number: AG20043	Calculated MW: 1615 aa, 179 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50882-4, LRP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60619-5-PBS. Detection antibody: 60619-3-PBS. Standard:Ag20043. Range: 0.391-100 ng/mL.