

For Research Use Only

ZBTB44 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60621-3-PBS



Basic Information

Catalog Number:

60621-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22369

GenBank Accession Number:

BC071729

GeneID (NCBI):

29068

UNIPROT ID:

Q8NCP5

Full Name:

zinc finger and BTB domain
containing 44

Calculated MW:

570 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B4G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50884-2, ZBTB44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60621-3-PBS. Detection antibody: 60621-2-PBS. Standard: Ag22369. Range: 0.781-100 ng/mL.