

For Research Use Only

# ZBTB44 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60621-5-PBS



## Basic Information

Catalog Number:

60621-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22369

GenBank Accession Number:

BC071729

GeneID (NCBI):

29068

UNIPROT ID:

Q8NCP5

Full Name:

zinc finger and BTB domain

containing 44

Calculated MW:

570 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C3C11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50884-3, ZBTB44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60621-4-PBS. Detection antibody: 60621-5-PBS. Standard:Ag22369. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50884-4, ZBTB44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60621-6-PBS. Detection antibody: 60621-5-PBS. Standard:Ag22369. Range: 0.781-100 ng/mL