

For Research Use Only

CYP19A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60624-3-PBS



Basic Information

Catalog Number: 60624-3-PBS	GenBank Accession Number: BC022896	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1588	CloneNo.: 1A6G4
Source: Mouse	UNIPROT ID: P11511	
Isotype: IgG1	Full Name: cytochrome P450, family 19, subfamily A, polypeptide 1	
Immunogen Catalog Number: AG9896	Calculated MW: 20 kDa, 25 kDa, 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50890-2, CYP19A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60624-3-PBS. Detection antibody: 60624-4-PBS. Standard: Ag9896. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50890-3, CYP19A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60624-3-PBS. Detection antibody: 60624-5-PBS. Standard: Ag9896. Range: 0.098-100 ng/mL