

For Research Use Only

SPATA17 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60625-1-PBS



Basic Information

Catalog Number: 60625-1-PBS	GenBank Accession Number: BC014608	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 128153	CloneNo.: 1H2E1
Source: Mouse	UNIPROT ID: Q96L03	
Isotype: IgG1	Full Name: spermatogenesis associated 17	
Immunogen Catalog Number: AG27909	Calculated MW: 361 aa, 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50892-1, SPATA17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60625-1-PBS. Detection antibody: 60625-2-PBS. Standard:Ag27909. Range: 3.125-200 ng/mL



Cytometric bead array standard curve of MP50892-2, SPATA17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60625-1-PBS. Detection antibody: 60625-3-PBS. Standard:Ag27909. Range: 3.125-200 ng/mL