

For Research Use Only

P16-INK4A Monoclonal antibody, PBS Only (Detector)

Catalog Number:60626-5-PBS



Basic Information

Catalog Number:

60626-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1328

GenBank Accession Number:

BC021998

GeneID (NCBI):

1029

UNIPROT ID:

Q8N726

Full Name:

cyclin-dependent kinase inhibitor 2A
(melanoma, p16, inhibits CDK4)

Calculated MW:

16 kDa

Purification Method:

Protein G purification

CloneNo.:

2F12E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51207-1, P16-INK4A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60626-4-PBS. Detection antibody: 60626-5-PBS. Standard:Ag1328. Range: 0.098-100 ng/mL