

For Research Use Only

BCDIN3D Monoclonal antibody, PBS Only (Detector)

Catalog Number:60628-2-PBS



Basic Information

Catalog Number: 60628-2-PBS	GenBank Accession Number: BC053560	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 144233	CloneNo.: 3G9A1
Source: Mouse	UNIPROT ID: Q7Z5W3	
Isotype: IgG1	Full Name: BCDIN3 domain containing	
Immunogen Catalog Number: AG17249	Calculated MW: 292 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50895-1, BCDIN3D Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60628-1-PBS. Detection antibody: 60628-2-PBS. Standard:Ag17249. Range: 0.391-100 ng/mL.