

For Research Use Only

# TMED9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60631-6-PBS



## Basic Information

Catalog Number:

60631-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13754

GenBank Accession Number:

BC001123

GeneID (NCBI):

54732

UNIPROT ID:

Q9BVK6

Full Name:

transmembrane emp24 protein  
transport domain containing 9

Calculated MW:

235 aa, 27 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2H3B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50898-4, TMED9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60631-1-PBS. Detection antibody: 60631-6-PBS. Standard: Ag13754. Range: 1.563-25 ng/mL