

For Research Use Only

MBTPS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60636-2-PBS



Basic Information

Catalog Number:

60636-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26396

GenBank Accession Number:

BC026330

GeneID (NCBI):

8720

UNIPROT ID:

Q14703

Full Name:

membrane-bound transcription factor
peptidase, site 1

Calculated MW:

1052 aa, 118 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A1B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50910-1, MBTPS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60636-1-PBS. Detection antibody: 60636-2-PBS. Standard:Ag26396. Range: 0.195-12.5 ng/mL



Cytometric bead array standard curve of MP50910-2, MBTPS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60636-3-PBS. Detection antibody: 60636-2-PBS. Standard:Ag26396. Range: 3.125-200 ng/mL