

For Research Use Only

ABCG1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60637-2-PBS



Basic Information

Catalog Number: 60637-2-PBS	GenBank Accession Number: BC029158	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9619	CloneNo.: 1E5A10
Source: Mouse	UNIPROT ID: P45844	
Isotype: IgG1	Full Name: ATP-binding cassette, sub-family G (WHITE), member 1	
Immunogen Catalog Number: AG17329	Calculated MW: 678 aa, 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50911-1, ABCG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60637-1-PBS. Detection antibody: 60637-2-PBS. Standard: Ag17329. Range: 0.391-100 ng/mL.