

For Research Use Only

ABCG1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60637-4-PBS



Basic Information

Catalog Number:

60637-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG17329

GenBank Accession Number:

BC029158

GeneID (NCBI):

9619

UNIPROT ID:

P45844

Full Name:

ATP-binding cassette, sub-family G
(WHITE), member 1

Calculated MW:

678 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

3B1A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50911-2, ABCG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60637-3-PBS. Detection antibody: 60637-4-PBS. Standard: Ag17329. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50911-3, ABCG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60637-5-PBS. Detection antibody: 60637-4-PBS. Standard: Ag17329. Range: 0.391-6.25 ng/mL.