

For Research Use Only

TSPAN4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60638-2-PBS



Basic Information

Catalog Number: 60638-2-PBS	GenBank Accession Number: BC000389	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 7106	CloneNo.: 1F8G6
Source: Mouse	UNIPROT ID: O14817	
Isotype: IgG1	Full Name: tetraspanin 4	
Immunogen Catalog Number: AG16938	Calculated MW: 238 aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50913-1, TSPAN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60638-1-PBS. Detection antibody: 60638-2-PBS. Standard: Ag16938. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50913-4, TSPAN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60638-5-PBS. Detection antibody: 60638-2-PBS. Standard: Ag16938. Range: 0.098-100 ng/mL