

For Research Use Only

SPAG5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60644-1-PBS



Basic Information

| | | |
|--|---|--|
| Catalog Number: 60644-1-PBS | GenBank Accession Number: BC000322 | Purification Method: Protein A Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 10615 | CloneNo.: 3D1G7 |
| Source: Mouse | UNIPROT ID: Q96R06 | |
| Isotype: IgG2b | Full Name: sperm associated antigen 5 | |
| Immunogen Catalog Number: AG6550 | Calculated MW: 134 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50922-1, SPAG5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60644-1-PBS. Detection antibody: 60644-2-PBS. Standard: Ag6550. Range: 6.25-200 ng/mL