

For Research Use Only

SPAG5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60644-3-PBS



Basic Information

Catalog Number: 60644-3-PBS	GenBank Accession Number: BC000322	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10615	CloneNo.: 1D10D10
Source: Mouse	UNIPROT ID: Q96R06	
Isotype: IgG2a	Full Name: sperm associated antigen 5	
Immunogen Catalog Number: AG6550	Calculated MW: 134 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50922-2, SPAG5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60644-3-PBS. Detection antibody: 60644-4-PBS. Standard: Ag6550. Range: 3.125-200 ng/mL.