

For Research Use Only

# AP1B1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60645-3-PBS



## Basic Information

Catalog Number:

60645-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10417

GenBank Accession Number:

BC046242

GeneID (NCBI):

162

UNIPROT ID:

Q10567

Full Name:

adaptor-related protein complex 1,  
beta 1 subunit

Calculated MW:

919 aa, 101 kDa

Purification Method:

Protein G purification

CloneNo.:

3A11A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50927-2, AP1B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60645-3-PBS. Detection antibody: 60645-4-PBS. Standard:Ag10417. Range: 3.125-100 ng/mL.