

For Research Use Only

Claudin 2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60650-1-PBS



Basic Information

Catalog Number:

60650-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1866

GenBank Accession Number:

BC014424

GeneID (NCBI):

9075

UNIPROT ID:

P57739

Full Name:

claudin 2

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F12C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50931-1, Claudin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60650-1-PBS. Detection antibody: 60650-2-PBS. Standard: Ag1866. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50931-3, Claudin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60650-1-PBS. Detection antibody: 60650-5-PBS. Standard: Ag1866. Range: 0.098-100 ng/mL