

For Research Use Only

FMO3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60651-3-PBS



Basic Information

Catalog Number: 60651-3-PBS	GenBank Accession Number: BC032016	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2328	CloneNo.: 1G12F3
Source: Mouse	UNIPROT ID: P31513	
Isotype: IgG1	Full Name: flavin containing monooxygenase 3	
Immunogen Catalog Number: AG11405	Calculated MW: 532 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50933-2, FMO3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60651-1-PBS. Detection antibody: 60651-3-PBS. Standard:Ag11405. Range: 0.098-100 ng/mL.