

For Research Use Only

TRIM36 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60653-2-PBS



Basic Information

Catalog Number:

60653-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG23227

GenBank Accession Number:

BC130334

GeneID (NCBI):

55521

UNIPROT ID:

Q9NQ86

Full Name:

tripartite motif-containing 36

Purification Method:

Protein A Magarose purification

CloneNo.:

2F5C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50936-1, TRIM36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60653-1-PBS. Detection antibody: 60653-2-PBS. Standard:Ag23227. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50936-2, TRIM36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60653-3-PBS. Detection antibody: 60653-2-PBS. Standard:Ag23227. Range: 0.391-50 ng/mL.