

For Research Use Only

ADAT2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60654-1-PBS



Basic Information

Catalog Number:

60654-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4717

GenBank Accession Number:

BC037955

GeneID (NCBI):

134637

UNIPROT ID:

Q7Z6V5

Full Name:

adenosine deaminase, tRNA-specific
2, TAD2 homolog (S. cerevisiae)

Calculated MW:

191 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H9G10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50937-1, ADAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60654-1-PBS. Detection antibody: 60654-2-PBS. Standard:Ag4717. Range: 1.563-100 ng/mL.