

For Research Use Only

# ADAT2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60654-4-PBS



## Basic Information

Catalog Number:

60654-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4717

GenBank Accession Number:

BC037955

GeneID (NCBI):

134637

UNIPROT ID:

Q7Z6V5

Full Name:

adenosine deaminase, tRNA-specific  
2, TAD2 homolog (S. cerevisiae)

Calculated MW:

191 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C4C7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50937-2, ADAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60654-3-PBS. Detection antibody: 60654-4-PBS. Standard:Ag4717. Range: 1.563-100 ng/mL.