

For Research Use Only

SPAM1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60655-2-PBS



Basic Information

Catalog Number:

60655-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1682

GenBank Accession Number:

BC026163

GeneID (NCBI):

6677

UNIPROT ID:

P38567

Full Name:

SPAM1

Calculated MW:

58 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G1E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50939-1, SPAM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60655-1-PBS. Detection antibody: 60655-2-PBS. Standard:Ag1682. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP50939-2, SPAM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60655-3-PBS. Detection antibody: 60655-2-PBS. Standard:Ag1682. Range: 3.125-100 ng/mL