

For Research Use Only

GPR161 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60657-2-PBS



Basic Information

Catalog Number: 60657-2-PBS	GenBank Accession Number: BC028163	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 23432	CloneNo.: 3C8B12
Source: Mouse	UNIPROT ID: Q8N6U8	
Isotype: IgG2b	Full Name: G protein-coupled receptor 161	
Immunogen Catalog Number: AG27121	Calculated MW: 529 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50941-1, GPR161 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60657-1-PBS. Detection antibody: 60657-2-PBS. Standard:Ag27121. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP50941-2, GPR161 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60657-3-PBS. Detection antibody: 60657-2-PBS. Standard:Ag27121. Range: 1.563-100 ng/mL