

For Research Use Only

# C3orf10 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60663-4-PBS



## Basic Information

<b>Catalog Number:</b> 60663-4-PBS	<b>GenBank Accession Number:</b> BC019303	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 55845	<b>CloneNo.:</b> 2E9D7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q8WUW1	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> chromosome 3 open reading frame 10	
<b>Immunogen Catalog Number:</b> AG33165	<b>Calculated MW:</b> 9 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50946-2, C3orf10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60663-3-PBS. Detection antibody: 60663-4-PBS. Standard: Ag33165. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50946-4, C3orf10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60663-4-PBS. Detection antibody: 60663-7-PBS. Standard: Ag33165. Range: 0.195-50 ng/mL