

For Research Use Only

DGAT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60664-1-PBS



Basic Information

Catalog Number: 60664-1-PBS	GenBank Accession Number: BC015762	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8694	CloneNo.: 1G2D1
Source: Mouse	UNIPROT ID: O75907	
Isotype: IgG2b	Full Name: diacylglycerol O-acyltransferase homolog 1 (mouse)	
Immunogen Catalog Number: AG32708	Calculated MW: 488 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50948-1, DGAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60664-1-PBS. Detection antibody: 60664-2-PBS. Standard:Ag32708. Range: 0.098-100 ng/mL.