

For Research Use Only

DGAT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60664-4-PBS



Basic Information

Catalog Number: 60664-4-PBS	GenBank Accession Number: BC015762	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 8694	CloneNo.: 2B6B4
Source: Mouse	UNIPROT ID: O75907	
Isotype: IgG1	Full Name: diacylglycerol O-acyltransferase homolog 1 (mouse)	
Immunogen Catalog Number: AG32708	Calculated MW: 488 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50948-2, DGAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60664-3-PBS. Detection antibody: 60664-4-PBS. Standard:Ag32708. Range: 0.195-100 ng/mL.