

For Research Use Only

DGAT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60664-6-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 60664-6-PBS | GenBank Accession Number: BC015762 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 8694 | CloneNo.: 2G5H12 |
| Source: Mouse | UNIPROT ID: O75907 | |
| Isotype: IgG1 | Full Name: diacylglycerol O-acyltransferase homolog 1 (mouse) | |
| Immunogen Catalog Number: AG32708 | Calculated MW: 488 aa, 55 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50948-3, DGAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60664-5-PBS. Detection antibody: 60664-6-PBS. Standard:Ag32708. Range: 0.781-100 ng/mL.