

For Research Use Only

# CORO6 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60665-3-PBS



## Basic Information

Catalog Number: 60665-3-PBS	GenBank Accession Number: BC064514	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 84940	CloneNo.: 1B10A2
Source: Mouse	UNIPROT ID: Q6QEF8	
Isotype: IgG1	Full Name: coronin 6	
Immunogen Catalog Number: AG11144	Calculated MW: 237aa,27 kDa; 286aa,32 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50950-2, CORO6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60665-3-PBS. Detection antibody: 60665-4-PBS. Standard:Ag11144. Range: 0.098-25 ng/mL