

For Research Use Only

IGF2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60668-2-PBS



Basic Information

Catalog Number:

60668-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21397

GenBank Accession Number:

BC000531

GeneID (NCBI):

3481

UNIPROT ID:

P01344

Full Name:

insulin-like growth factor 2
(somatomedin A)

Calculated MW:

236 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H9B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50953-1, IGF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60668-1-PBS. Detection antibody: 60668-2-PBS. Standard: Ag21397. Range: 0.391-200 ng/mL.