

For Research Use Only

PRPF4B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60669-2-PBS



Basic Information

Catalog Number:

60669-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG11006

GenBank Accession Number:

BC034969

GeneID (NCBI):

8899

UNIPROT ID:

Q13523

Full Name:

PRP4 pre-mRNA processing factor 4
homolog B (yeast)

Calculated MW:

1007 aa, 117 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D1E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50956-1, PRPF4B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60669-1-PBS. Detection antibody: 60669-2-PBS. Standard:Ag11006. Range: 0.098-100 ng/mL.