

For Research Use Only

IL18BP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60671-4-PBS



Basic Information

Catalog Number: 60671-4-PBS	GenBank Accession Number: BC044215	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10068	CloneNo.: 1F3B3
Source: Mouse	UNIPROT ID: O95998	
Isotype: IgG2b	Full Name: interleukin 18 binding protein	
Immunogen Catalog Number: AG5638	Calculated MW: 194 aa, 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50961-3, IL18BP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60671-4-PBS. Detection antibody: 60671-5-PBS. Standard:Ag5638. Range: 0.781-100 ng/mL.