

For Research Use Only

NUS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60672-2-PBS



Basic Information

Catalog Number:

60672-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8959

GenBank Accession Number:

BC013026

GeneID (NCBI):

116150

UNIPROT ID:

Q96E22

Full Name:

nuclear undecaprenyl pyrophosphate
synthase 1 homolog (S. cerevisiae)

Calculated MW:

293 aa, 33 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C1F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50965-1, NUS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60672-1-PBS. Detection antibody: 60672-2-PBS. Standard:Ag8959. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50965-2, NUS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60672-3-PBS. Detection antibody: 60672-2-PBS. Standard:Ag8959. Range: 0.098-100 ng/mL.