

For Research Use Only

NUS1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60672-4-PBS



Basic Information

Catalog Number: 60672-4-PBS	GenBank Accession Number: BC013026	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 116150	CloneNo.: 2A7B10
Source: Mouse	UNIPROT ID: Q96E22	
Isotype: IgG1	Full Name: nuclear undecaprenyl pyrophosphate synthase 1 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG8959	Calculated MW: 293 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50965-3, NUS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60672-4-PBS. Detection antibody: 60672-5-PBS. Standard: Ag8959. Range: 0.098-100 ng/mL.