

For Research Use Only

# ZMPSTE24 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60675-1-PBS



## Basic Information

Catalog Number:

60675-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG33613

GenBank Accession Number:

BC037283

GeneID (NCBI):

10269

UNIPROT ID:

O75844

Full Name:

zinc metallopeptidase (STE24 homolog, *S. cerevisiae*)

Calculated MW:

475 aa, 55 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B5C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50967-1, ZMPSTE24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60674-1-PBS. Detection antibody: 60675-1-PBS. Standard: Ag33613. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50967-2, ZMPSTE24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60674-2-PBS. Detection antibody: 60675-1-PBS. Standard: Ag33613. Range: 0.098-12.5 ng/mL