

For Research Use Only

SREBF2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60677-3-PBS



Basic Information

Catalog Number: 60677-3-PBS	GenBank Accession Number: BC056158	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6721	CloneNo.: 1A3A5
Source: Mouse	UNIPROT ID: Q12772	
Isotype: IgG1	Full Name: sterol regulatory element binding transcription factor 2	
Immunogen Catalog Number: AG28241	Calculated MW: 124 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50973-2, SREBF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60677-3-PBS. Detection antibody: 60677-4-PBS. Standard:Ag28241. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50973-3, SREBF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60677-3-PBS. Detection antibody: 60677-5-PBS. Standard:Ag28241. Range: 0.098-100 ng/mL.