

For Research Use Only

# SREBF2 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60677-5-PBS



## Basic Information

<b>Catalog Number:</b> 60677-5-PBS	<b>GenBank Accession Number:</b> BC056158	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6721	<b>CloneNo.:</b> 1B6F7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q12772	
<b>Isotype:</b> IgG1	<b>Full Name:</b> sterol regulatory element binding transcription factor 2	
<b>Immunogen Catalog Number:</b> AG28241	<b>Calculated MW:</b> 124 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50973-3, SREBF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60677-3-PBS. Detection antibody: 60677-5-PBS. Standard:Ag28241. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50973-4, SREBF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60677-2-PBS. Detection antibody: 60677-5-PBS. Standard:Ag28241. Range: 0.195-100 ng/mL.