

For Research Use Only

PSEN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60678-3-PBS



Basic Information

Catalog Number: 60678-3-PBS	GenBank Accession Number: BC011729	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5663	CloneNo.: 3E11F5
Source: Mouse	UNIPROT ID: P49768	
Isotype: IgG1	Full Name: presenilin 1	
Immunogen Catalog Number: AG14464	Calculated MW: 467 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50976-2, PSEN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60678-3-PBS. Detection antibody: 60678-4-PBS. Standard: Ag14464. Range: 0.391-100 ng/mL.