

For Research Use Only

# SPAG8 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60679-2-PBS



## Basic Information

Catalog Number:

60679-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5150

GenBank Accession Number:

BC019247

GeneID (NCBI):

26206

UNIPROT ID:

Q99932

Full Name:

sperm associated antigen 8

Calculated MW:

45 kDa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C6H5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50977-1, SPAG8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60679-1-PBS. Detection antibody: 60679-2-PBS. Standard: Ag5150. Range: 1.563-100 ng/mL.